

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курган (3522)50-90-47
Курск (4712)77-13-04
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Ноябрьск (3496)41-32-12

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Пермь (342)205-81-47
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саранск (8342)22-96-24
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70


Казахстан +7(7172)727-132


Киргизия +996(312)96-26-47

mvi@nt-rt.ru || <https://mgb.nt-rt.ru>


Троакары с шарико-магнитными, силиконовыми, трубчатыми, откидными клапанами, троакары большого диаметра, безопасные троакары Техническое описание


MAGNETIC BALL – TROCAR SLEEVES



Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-24055	5.5 mm	95.0 mm	Trocar Sleeve with Magnetic Ball ■ Stainless steel
	429-24100	10.0 mm	97.0 mm	Trocar Sleeve with Magnetic Ball ■ Stainless steel

SILICONE VALVE TROCAR-SLEEVE



Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-25055	5.5 mm	95.0 mm	Trocar Sleeve with Silicone Valve ■ Stainless steel

TRUMPET VALVE – TROCAR SLEEVE

Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-40005	5.5 mm	95.0 mm	Metal Shaft ■ Straight
	429-40013	10.0 mm	97.0 mm	
	429-40014	11.0 mm	101.0 mm	

FLAP VALVE – TROCAR SLEEVES

Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-40201	3.5 mm	70.0 mm	Metal Shaft ■ Straight
	429-45201	5.0 mm	92.0 mm	
	429-40001	5.5 mm	95.0 mm	
	429-40004	7.0 mm	95.0 mm	
	429-45001	10.0 mm	97.0 mm	
	429-40002	11.0 mm	101.0 mm	
	429-40003	12.5 mm	98.0 mm	
	429-40011	5.5 mm	95.0 mm	Metal Shaft ■ With retaining thread
	429-40012	11.0 mm	101.0 mm	
	429-41012	12.5 mm	98.0 mm	
	429-40051	5.5 mm	95.0 mm	Plastic Shaft ■ Straight
	429-40052	11.0 mm	101.0 mm	
	429-40054	12.5 mm	98.0 mm	
	429-40061	5.5 mm	95.0 mm	Plastic Shaft ■ With retaining thread
	429-40062	11.0 mm	101.0 mm	
	429-40064	12.5 mm	98.0 mm	




FLAP VALVE – TROCAR SLEEVES

Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-40202	3.5 mm	70.0 mm	Metal Shaft ■ Straight ■ With manual opening
	429-40020	5.0 mm	92.0 mm	
	429-40021	5.5 mm	95.0 mm	
	429-40025	10.0 mm	97.0 mm	
	429-40022	11.0 mm	101.0 mm	
	429-40024	12.5 mm	98.0 mm	
	429-40021/15	5.5 mm	150.0 mm	Bariatric Surgery
	429-40022/15	11.0 mm	150.0 mm	
	429-40024/15	12.5 mm	150.0 mm	
	429-40031	5.5 mm	95.0 mm	Metal Shaft ■ With retaining thread ■ With manual opening
	429-40032	11.0 mm	101.0 mm	
	429-40034	12.5 mm	98.0 mm	
	429-40071	5.5 mm	95.0 mm	Plastic Shaft ■ Straight ■ With manual opening
	429-40072	11.0 mm	101.0 mm	
	429-40074	12.5 mm	98.0 mm	
	429-40081	5.5 mm	95.0 mm	Plastic Shaft ■ With retaining thread ■ With manual opening
	429-40082	11.0 mm	101.0 mm	
	429-40084	12.5 mm	98.0 mm	




SILICONE VALVE – TROCAR SLEEVES

Products	With Luer Stopcock	Diameter	Working Length	Description
	Order-No.			
	429-40280	3.5 mm	70.0 mm	Metal Shaft ■ Straight
	429-40281	5.0 mm	92.0 mm	
	429-40282	5.5 mm	95.0 mm	
	429-40283	10.0 mm	97.0 mm	
	429-40284	11.0 mm	101.0 mm	
	429-40285	12.5 mm	98.0 mm	
	429-40306	5.5 mm	95.0 mm	Metal Shaft ■ With retaining thread
	429-40308	11.0 mm	101.0 mm	
	429-40309	12.5 mm	98.0 mm	
	429-40311	5.5 mm	95.0 mm	Plastic Shaft ■ Straight
	429-40313	11.0 mm	101.0 mm	
	429-40314	12.5 mm	98.0 mm	
	429-40321	5.5 mm	95.0 mm	Plastic Shaft ■ With retaining thread
	429-40323	11.0 mm	101.0 mm	
	429-40324	12.5 mm	98.0 mm	





BIG DIAMETER FLAP VALVE – TROCAR SLEEVES WITH MANUAL OPENING

Products	With	Diameter	Working Length	Description
	Luer Stopcock			
	Order-No.			
	429-40041	15.0 mm	110.0 mm	Metal Shaft
	429-40042	20.0 mm	120.0 mm	■ Straight
	429-23900			Abdominal Speculum

TROCARS

Products	Order-No.	Diameter	For Sleeve Working Length	Description
	429-40203	3.5 mm	70.0 mm	Trocar
	429-40130	5.0 mm	92.0 mm	■ Conical
	429-40131	5.5 mm	95.0 mm	
	429-40170	7.0 mm	95.0 mm	
	429-40134	10.0 mm	97.0 mm	
	429-40133	11.0 mm	101.0 mm	
	429-40137	12.5 mm	98.0 mm	
	429-40131/15	5.5 mm	150.0 mm	
	429-40133/15	11.0 mm	150.0 mm	Bariatric Surgery
	429-40137/15	12.5 mm	150.0 mm	
	429-40204	3.5 mm	70.0 mm	Trocar
	429-40136	5.0 mm	92.0 mm	■ Pyramidal
	429-40141	5.5 mm	95.0 mm	
	429-40159	7.0 mm	95.0 mm	
	429-40138	10.0 mm	97.0 mm	
	429-40143	11.0 mm	101.0 mm	
	429-40139	12.5 mm	98.0 mm	
	429-40141/15	5.5 mm	150.0 mm	
	429-40143/15	11.0 mm	150.0 mm	Bariatric Surgery
	429-40139/15	12.5 mm	150.0 mm	
	429-40146	3.5 mm	70.0 mm	Trocar
	429-40147	5.0 mm	92.0 mm	■ Conical blunt
	429-40148	5.5 mm	95.0 mm	
	429-40160	7.0 mm	95.0 mm	
	429-40149	10.0 mm	97.0 mm	
	429-40151	11.0 mm	101.0 mm	
	429-40152	12.5 mm	98.0 mm	
	429-40153	5.0 mm	92.0 mm	Trocar
	429-40154	5.5 mm	95.0 mm	■ Pyramidal with hollow grinding
	429-40163	7.0 mm	95.0 mm	
	429-40135	10.0 mm	97.0 mm	
	429-40156	11.0 mm	101.0 mm	
	429-40155	12.5 mm	98.0 mm	

SAFETY TROCARS

Products	Order-No.	Diameter	For Sleeve Working Length	Description
	Order-No.			
	429-40103	5.0 mm	92.0 mm	Trocar ■ Pyramidal
	429-40101	5.5 mm	95.0 mm	
	429-40104	10.0 mm	97.0 mm	
	429-40102	11.0 mm	101.0 mm	
	429-40105	12.5 mm	98.0 mm	
	429-40113	5.0 mm	92.0 mm	Trocar ■ Half-moon
	429-40111	5.5 mm	95.0 mm	
	429-40114	10.0 mm	97.0 mm	
	429-40112	11.0 mm	101.0 mm	
	429-40115	12.5 mm	98.0 mm	
	429-40125	5.5 mm	95.0 mm	Trocar ■ V-style
	429-40126	10.0 mm	97.0 mm	
	429-40127	11.0 mm	101.0 mm	
	670-95341 (100 mm)			Verres Needle Ø 2 mm
	670-95342 (120 mm)			
	670-95340 (150 mm)			
	670-95351 (100 mm)			Verres Needle Ø 2 mm
	670-95352 (120 mm)			
	670-95350 (150 mm)			



По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курган (3522)50-90-47
Курск (4712)77-13-04
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Ноябрьск (3496)41-32-12

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Пермь (342)205-81-47
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саранск (8342)22-96-24
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

mvi@nt-rt.ru || <https://mgb.nt-rt.ru>